- EN: Instructions for Use**
Fustar™ Steerable Introducer
- FR: Mode d'emploi**
Fustar™ Introducteur orientable
- DE: Gebrauchsanweisung**
Fustar™ Steuerbares Einführungsinstrument
- ES: Instrucciones de uso**
Fustar™ Introducutor dirigitible
- PT: Instruções de utilização**
Fustar™ Introduutor orientável
- TR: Kullanım Talimatları**
Fustar™ Yönlendirilebilir Uygulayıcı
- IT: Istruzioni per l'uso**
Introduttore orientabile Fustar™
- DA: Brugsanvisning**
Fustar™ Justerbar buet kappe
- FI: Käyttöohjeet**
Fustar™ Säädettävä kaareva vaippa
- NL: Gebruiksaanwijzing**
Fustar™ Justerbar buet kappe
- SV: Bruksanvisning**
Fustar™ Justerbar krökt mantel
- EL: Οδηγίες χρήσης**
Κατευθυνόμενος εισαγωγέας Fustar™
- RU: Инструкция по применению**
Управляемый интродьюсер Fustar™


0344

	Catalogue number Numéro de catalogue Katalognummer Número de catálogo Número de catálogo Katalog numarası Numero di catalogo Katalognummer Luettelonumero Catalogusnummer Katalognummer Αριθμόςκαταλόγου Номер по каталогу
	Use by À utiliser jusqu'au Zu verwenden bis einschließlich No utilizar después de Não utilizar depois de Son kullanım tarihi Utilizzare entro Kan anvendes til og med Käytettävä viimeistään Uiterste gebruiksdatum Får användas till och med Χρήσήςως Использовать до
	Date of manufacture Date de fabrication Herstellungsdatum Fecha de fabricación Data de fabrico Üretim tarihi Data di fabbricazione Fabrikationsdato Valmistuspäivämäärä Productiedatum Tillverkningsdatum Ημερομηνίακατασκευής Дата изготовления
	Serial Number Numéro de série Seriennummer Número de serie Número de série Seri numarası Numero di serie Seriennummer Sarjnumero Seriennummer Seriennummer Αριθμόςσειράς Серийный номер
	Lot Number Numéro de lot Chargennummer Número de lote Número de lote Lot Numarası Numero di lotto Partinummer Eränumero Partijnummer Lotnummer Αριθμός παρτίδας Номер партии
	Consult instructions for use Consulter le mode d'emploi Gebrauchsanweisung beachten Consultar las instrucciones de uso Consultar instruções de utilização Kullanım talimatına bakın Consultare le istruzioni per l'uso Se brugsanvisningen Katso käyttöohjeet Raadpleeg gebruiksaanwijzing Se bruksanvisning Συμβουλευτείτε τις Οδηγίες χρήσης Обратитесь к инструкции по применению

	Sterilized using ethylene oxide Stérilisation par oxyde d'éthylène Sterilisiert mit ethylenoxid Esterilizado mediante óxido de etileno Etilen oksit ile sterilize edilmiştir Sterilizzato a ossido di etilene Steriliseret med etylenoxid Steriloitu etyleenioksidilla Gesteriliseerd met ethyleenoxide Steriliserad med etylenoxid Αποστειρωμένο με αιθυλενοξείδιο Стерилизовано этиленоксидом
	Non-pyrogenic Apyrogène Pyrogenfrei Apirógeno Não pirogênico Pirojenik değildir Apirogeno Ikke-pyrogen Pyrogeeniton Niet-pyrogeen Pyrogenfri Μη πυρετογόνο Аπυρογενно
	Do not reuse Ne pas réutiliser Nicht wiederverwenden No reutilizar Não reutilizável Yeniden kullanmayın Non riutilizzare Må ikke genbruges Älä käytä uudelleen Niet opnieuw gebruiken Får ej återanvändas Μη επαναχρησιμοποιείτε Не использовать повторно
	Do not use if package is damaged Ne pas utiliser si l'emballage est endommagé Nicht verwenden, wenn die Verpackung beschädigt ist No utilizar si el envase está dañado Não use se a embalagem estiver danificada Ambalaj hasar gördüyse kullanmayın Non utilizzare se l'imballo non è integro NMå ikke anvendes, hvis emballagen er beskadiget Älä käyttää, jos pakkaus on vaurioitunut Niet gebruiken als de verpakking beschadigd is Får ej användas om förpackningen är skadad Μη χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά Не использовать в случае повреждения упаковки
	Do not resterilize Ne pas restériliser Nicht reesterilisieren No reesterilizar Não reesterilizar Yeniden sterilize etmeyin Non risterilizzare Må ikke reesteriliseres Älä steriloitua uudelleen Niet opnieuw steriliseren Får ej omsteriliseras Μη επαναποστειρώνετε Не стерилизовать повторно

	Store at room temperature Conserver à température ambiante Bei Raumtemperatur lagern Almacenar el producto a temperatura ambiente Guardar à temperatura ambiente Oda sıcaklığında saklayın Conservare a temperatura ambiente Skal opbevares ved stuetemperatur Säilytä huoneenlämmössä Bij kamertemperatuur bewaren Förvaras i rumstemperatur. Φυλάσσετε σε θερμοκρασία δωματίου Хранить при комнатной температуре
	Keep dry À conserver dans un endroit sec Trocken aufbewahren Manterer seco Manterha seco Kuru yerde saklayın Manterere asciutto Skal opbevares tørt Säilytä kuivassa Droog bewaren Förvaras tort Διατηρήστε στεγνό Хранить в сухом месте
	Keep away from sunlight Ne pas exposer à la lumière du soleil Vor Sonnenlicht schützen Manterer alejado de la luz del sol Manter afastado da luz do sol Güneş ışığından uzak tutun Non esporre ai raggi solari Må ikke udsættes for sollys Suojaa auringonvalolta Uit de buurt van zonlicht houden Får inte utsättas för solljus Διατηρείτε μακριά από το ηλιακό φως Беречь от солнечных лучей
	Manufacturer Fabricant Hersteller Fabricante Fabricante İmalatçı Fabbicante Fabrikant Valmistaja Fabrikant Tillverkare Κατασκευαστής Производитель

	Authorized representative in the European Community Représentant autorisé dans la communauté européenne Autorisierter Repräsentant in der Europäischen Gemeinschaft Representante autorizado en la Comunidad Europea Representante autorizado na Comunidade Europeia Avrupa Topluluğu'ndaki yetkili temsilci Rappresentante autorizzato nella Comunità europea Autoriseret repræsentant i EU Valtuutettu edustaja Euroopan yhteisön alueella Geautoriseerd vertegenwoordiger in de Europese Gemeenschap Auktoriserad representant inom EG Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα Уполномоченный представитель в Европейском сообществе
	Medical device Dispositif médical Medizinisches Gerät Producto sanitario Dispositivo médico Tıbbi Cihaz Dispositivo medico Medicinsk udstyr Lääkinnällinen laite Medisch apparaat Medicinteknisk produkt Ιατροτεχνολογικό προϊόν Медицинское оборудование
	Importer Importateur Importeur Importador Importador İthalatçı Importatore Importør Maahantuoja Importeur Importör Εισαγωγέας Импортёр
	Made in china Fabriqué en Chine Hergestellt in China Hecho en China Fabricado na China Çin'de üretilmiştir Fabricato in Cina Fremstillet i Kina Valmistettu Kiinassa Gemaakt in China Tillverkad i Kina Κατασκευασμένο στην Κίνα Сделано в КНР
	Unique device identifier Identifiant unique du dispositif Eindeutige Gerätekennung Identificación única del producto Identificador único do dispositivo (UDI) Benzersiz cihaz tanımlayıcısı Identificativo univoco del dispositivo unik udstyrsidentifikation Yksilöllinen laitetunniste Uniek apparaatnummer Unik produktidentifierring Αποκλειστική ταυτοποίηση ιατροτεχνολογικού προϊόντος Уникальный идентификатор устройства

Double sterile barrier system | Système de double
barrière stérile | System mit doppelter Sterilbarriere |
Sistema de doble barrera estéril | Sistema de barreira
estéril dupla | Çifte steril bariyer sistemi | Sistema a
doppia barriera sterile | System med steril dobbelt
barriere | Kaksinkertainen steriili estejärjestelmä |
Systeem met dubbele, steriele barrière | Dubbelt sterilt
barriärsystem | Σύστημα φραγμού διπλής αποστείρωσης
| Двойная стерильная барьерная система

1 DEVICE DESCRIPTION

The Fustar™ Steerable introducer comprises:

- 1) Deflectable sheath (Fig.1 & Fig.2)
- 2) Dilator (Fig.3 & Fig.4)
- 3) Haemostatic valve (Fig.5)
- 4) Loader (Fig.6)
- 5) Delivery cable (optional) (Fig.7)

1.1 Sheath

- Radiopaque sheath and a radiopaque marker band incorporated on the distal end of the sheath.
- The distal end of the sheath could be deflected by rolling the knob located on the handle. When rolling the knob clockwise, the distal end of sheath will deflect to larger angle. Oppositely, when rolling the knob counter clockwise, the distal end of sheath will recover its original status.
- The product is divided into two series and the details are shown in Section 2.


Fig.1 Steerable introducer(800 series)


Fig.2 Steerable introducer(550/700/900 series)

1.2 Dilator

- The dilator is used to ease penetration of tissues and the vessel wall.


Fig.3 Dilator(800 series)


Fig.4 Dilator(550/700/900 series)

1.3 Haemostatic valve

- The haemostatic valve is used to connect with the hub of the loader to prevent the blood leakage and to flush devices. Before flushing the sheath, shut down the stopcock clockwise to prevent ingress of air. When the device enters the sheath, open the stopcock counter-clockwise to push the delivery cable.


Fig.5 Haemostatic valve

1.4 Loader

- The loader is used to load an occlusion device(if applicable) into the sheath and connect it to the sheath lumen.


Fig.6 Loader

1.5 Delivery cable (optional)

- The delivery cable is used to attach and deliver the Lifetech Occluder, Vena Cava Filter, or Vascular Plug.


Fig.7 Delivery cable

2 PRODUCT SPECIFICATION

2.1 Sheath types

The sheath has various lumen diameters, deflectable section lengths and effective lengths of sheath as follows:

- There are 8 sizes of the lumen diameters ranging from 5F to 14F.
- The length of the deflectable distal end is either 30mm or 50mm.
- The effective lengths of the sheath are 550mm, 700mm, 800mm and 900mm.
- The Fustar™ Steerable Introducer is divided into two series according to the effective lengths, i.e. 550/700/900 series and 800 series.

2.2 550/700/900 series

This series contains a sheath and a dilator only. The details are shown in table.1.

Table.1 The product family of Fustar™ Steerable Introducer (550/700/900)

Curve Type	Specification			ID/mm	OD/mm	Deflectable Angle/°
	550	700	900			
Type A (30mm)	SVA5F-550	SVA5F-700	SVA5F-900	1.8	2.6	0~180
	SVA6F-550	SVA6F-700	SVA6F-900	2.1	2.9	
	SVA7F-550	SVA7F-700	SVA7F-900	2.5	3.3	
	SVA8F-550	SVA8F-700	SVA8F-900	2.8	3.6	
	SVA9F-550	SVA9F-700	SVA9F-900	3.1	3.9	
	SVA10F-550	SVA10F-700	SVA10F-900	3.5	4.3	
Type B (50mm)	SVB5F-550	SVB5F-700	SVB5F-900	1.8	2.6	0~180
	SVB6F-550	SVB6F-700	SVB6F-900	2.1	2.9	
	SVB7F-550	SVB7F-700	SVB7F-900	2.5	3.3	
	SVB8F-550	SVB8F-700	SVB8F-900	2.8	3.6	
	SVB9F-550	SVB9F-700	SVB9F-900	3.1	3.9	
	SVB10F-550	SVB10F-700	SVB10F-900	3.5	4.3	
	SVB12F-550	SVB12F-700	SVB12F-900	4.1	5.0	
						0~90

Curve Type	Specification			ID/mm	OD/mm	Deflectable Angle/°
	550	700	900			
	SVB14F-550	SVB14F-700	SVB14F-900	4.8	5.7	

Note: 1F=1 French=1/3 mm

2.3 800 series

This series contains a sheath, a dilator, a loader, a haemostatic valve and a delivery cable. The details are shown in table.2.

Table.2 The product family of Fustar™ Steerable Introducer(800)

Curve Type	Specification	ID/mm	OD/mm	Deflectable Angle/°
	800			
Type A (30mm)	SVA5F-800	1.8	2.6	0~180
	SVA6F-800	2.1	2.9	
	SVA7F-800	2.5	3.3	
	SVA8F-800	2.8	3.6	
	SVA9F-800	3.1	3.9	
	SVA10F-800	3.5	4.3	
Type B (50mm)	SVB12F-800	4.1	5.0	0~90
	SVB14F-800	4.8	5.7	

Note: 1F=1 French=1/3 mm

3 INDICATIONS AND USAGE

The Fustar™ Steerable Introducer is designed for the percutaneous interventional vascular procedure. The Fustar™ Steerable Introducer would provide access for the introduction and retrieval of the interventional device, such as occluder, balloon, vena cava filter, stent, diagnostic and guiding catheters or other device, fluid or media for diagnosis or intervention.

4 CONTRAINDICATIONS

Contraindications would be similar to those in which a standard introducer or delivery sheath could be used.

5 PRECAUTIONS

- 1) The Fustar™ Steerable Introducer is intended for use only by physicians trained and experienced with diagnostic and interventional techniques.
- 2) Do not attempt to use the Fustar™ Steerable Introducer before reading this instruction booklet and understanding the device and its use.
- 3) The product has been sterilized by EO. Do not use if the sterile barrier has been compromised in any way. Please carefully inspect all packages before use.

- 4) Before removing or inserting device through the introducer, aspirate through the side-arm of the valve to clear the introducer, and then flush with heparinized saline.
- 5) Remove and replace the introducer if the distal end of sheath could not be deflected.
- 6) The maximum diameter of the instrument or catheter to be introduced should be determined to ensure its safe passage through the introducer. All devices or catheters used with this product should move freely through the valve and sheath. Damage to the valve of introducer may result if the fit is tight.
- 7) When inserting, manipulating or withdrawing a device through an introducer, always maintain introducer position.
- 8) Do not attempt to insert or withdraw the guide wire and/ or introducer if resistance is felt.
- 9) Do not release the device from the delivery cable if the device does not conform to its original configuration or if the device position is unstable. Recapture the device and redeploy. If still unsatisfactory, recapture the device and replace with a new device.
- 10) Withdraw the dilator 80 mm back from the tip, when steering the curve of the distal tip for directional change. Take caution not to withdraw the dilator too quickly, or air embolism may occur.
- 11) For disposable, single-use only. Do not reuse, reprocess, or re-sterilize.
- 12) Disposal of product and packaging should be in accordance with hospital, administrative and/or local government policy.

6 POTENTIAL ADVERSE EVENTS

Possible adverse events may occur during the use of this device and include, but are not limited to:

- 1) Air embolism
- 2) Infection
- 3) Intimal tear
- 4) Hematoma
- 5) Perforation
- 6) Thrombus formation

7 DIRECTIONS FOR USE

Materials recommended for use with the delivery system

exchange-length 0.038-inch guidewire

Sheath Introduction

- 1) Upon removal from package, ensure the inner diameter(ID) of the sheath is appropriate for the maximum diameter of the device.
- 2) Using the side-arm of the valve, flush the steerable sheath by filling the steerable

sheath completely with heparinized saline. Meanwhile, flush the dilator with heparinized saline.

- 3) Insert the dilator completely into the steerable sheath, turn the rotating luer on the dilator clockwise to lock the components together, and then deliver both into the body via the puncture point along the guidewire.
- 4) When the sheath/dilator combination arrives at the target area, withdraw the dilator approximately 80 mm from the tip for easy deflection. Rotate the handle in a clockwise manner to bend the end of sheath. Manipulate the sheath tip to the target area.
- 5) Take out the guidewire and dilator slowly. Aspirate and flush the side-arm haemostatic valve.
- 6) Insert proper device as needed.

Sheath Removal

- 1) Insert the guidewire at least 10cm past the tip of the sheath.
- 2) Insert the dilator over the wire into the sheath as needed.
- 3) Rotate the knob counter-clockwise to recover the angle at the end of the sheath. Then withdraw the sheath from the target location to prevent from injuring the vascular wall.
- 4) Remove the guide wire.

8 EXPIRY DATE

Please use the device before the "use by" date specified on the label. The sterilization date is defined as the date of manufacture.

9 STORAGE AND TRANSPORTATION

- 1) The Fustar™ Steerable Introducer shall be kept in protective, well ventilated and clean room with no aggressive gas. The maximum limit for relative humidity is 80%.
- 2) Be careful in its handling during transportation. Do not allow the expose the sheath to wetness, unusual heat, or pressure during transportation.
- 3) For any question before using the Fustar™ Steerable Introducer, contact the sales agent or the manufacturer directly.

10 PACKAGE AND LABEL

The Fustar™ Steerable introducer is sterile by EO. All the components are sealed in two Tyvek pouches on which a label and a sterilization indicator are placed, then immobilized and protected on a scale board. The product is sterilized and put in a box with IFU, patient card, customer feedback form and certificate of compliance.

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Lifetech Scientific (Shenzhen) Co., Ltd