SeQure™

Snare System	EN 1-8
نظام السنار تعليمات الاستخدام	AR9-15
Ендоваскуларна примка за извличане на чужди тела Инструкции за употреба	BG16-24
Smyčkový system Návod k použití	CS 25-32
Snaresystem Brugsanvisning	DA 33-40
Schlingenkatheter-System Gebrauchsanweisung	DE41-49
Σύστημα αφαίρεσης με βρόχο Οδηγίες χρήσης	EL50-58
Sistema de extracción Instrucciones de uso	ES59-66
Võõrkeha haaramise ling Kasutusjuhend	ET 67-74
Silmukkajärjestelmä Käyttöohje	FI 75-82
Système d'extraction à l'anse Mode d'emploi	FR83-91
Sustav sa zamkom Upute za uporabu	HR92-99

Snare System Használati utasítás	HU 100-107
Sistema ad ansa Istruzioni per l'uso	IT 108-115
„Snare System“ Naudojimo instrukcijos	LT 116-123
Snaarsysteem Gebruiksaanwijzing	NL 124-131
Slyngekatetersystem Bruksanvisning	NO132-139
System końcówek do szycia Instrukcja użytkowania	PL 140-148
Sistema de captura (laço) Instruções de utilização	PT149-156
Sistem cateter cu laț Instrucțiuni de utilizare	RO 157-164
Система петлевой ловушки Инструкция по применению	RU 165-173
Slučkový system Návod na použitie	SK174-181
Sistem katetra z zanko Navodila za uporabo	SL182-189
Snare sistem Uputstvo za upotrebu	SR 190-197
Slyngsystem Bruksanvisning	SV 198-205
Kement Sistemi Kullanım Talimatları	TR206-213

	Date of manufacture تاريخ التصنيع Дата на производство Datum výroby Produktionsdato Herstellungsdatum Ημερομηνία κατασκευής Fecha de fabricación Tootmiskuu päev Valmistuspäivämäärä Date de fabrication Datum proizvodnje Gyártás dátuma Data di fabbricazione Pagaminimo data Datum van vervaardiging Produksjonsdato Data produkcji Data de fabric Data fabricației Дата изготовления Datum výroby Datum proizvodnje Datum proizvodnje Tillverkningsdatum Üretim tarihi
	Manufacturer الشركة المصنعة Производител Výrobce Producent Hersteller Κατασκευαστής Fabricante Tootja Valmistaja Fabricant Proizvođač Gyártó Fabbicante Gamintojas Fabrikant Produsent Producent Fabricante Producător Производитель Výrobca Proizvajalec Proizvođač Tillverkare Üretici
	Use-by date تاريخ انتهاء الصلاحية Срок на годност Datum spotřeby Sidste anvendelsesdato Verfallsdatum Ημερομηνία λήξης Fecha de caducidad Kõlblikusaeg Viimeinen käyttöpäivä Date de péremption Rok upotrebe Felhasználhatósági idő Data di scadenza Tinka naudoti iki Houdbaarheidsdatum Utløpsdato Termin ważności Prazo de validade A se utiliza înainte de Дата истечения срока годности Použite do Uporabite do datuma Rok upotrebe Förbrukningsdatum Son kullanma tarihi
	Serial number الرقم التسلسلي Сериен номер Sériové číslo Serienummer Seriennummer Αριθμός σειράς Número de serie Seerianumber Sarjanumero Numéro de série Serijski broj Sorozatszám Numero di serie Serijos numeris Serienummer Serienummer Numer seryjny Número de série Număr de serie Сериальный номер Sériové číslo Serijska številka Serijski broj Serienummer Seri numarası

LOT	Batch code كود التشغيلة Код на партида Kód šarže Batch-kode Chargennummer Αριθμός παρτίδας Código de lote Partiikood Eräkoodi Code de lot Šifra lota Sarzs-kód Codice del lotto Partijos kodas Batchcode Batchkode Kod partii Código do lote Cod lot Код партии Kód šarže Koda serije Kod grupe Satskod Parti kodu
STERILE EO	Sterilized using Ethylene Oxide معقم باستخدام أكسيد الإيثيلين Стерилизовано с этилен оксид Sterilizováno etylenoxidem Steriliseret med ethylenoxid Mit Ethylenoxid sterilisiert Αποστειρωμένο με τη χρήση οξειδίου του αιθυλενίου Esterilizado con óxido de etileno Steriliseeritud etüleenoksiidiga Steriloitu eteenioksidilla Stérilisé à l'oxyde d'éthylène Sterilizirano etilen-oksidom Etílenoxiddal sterilizálva Sterilizzato con ossido di etilene Sterilizuota etileno oksidu Gesteriliseerd met ethyleenoxide Sterilisert med etylenoksid Wyjałowiono tlenkiem etylenu Esterilizado utilizando óxido de etileno Sterilizat cu oxid de etilenă Стерилизовано этиленоксидом Sterilizované etylénoxidom Sterilizirano z etilen oksidom Sterilizovano pomoću etilen oksida Steriliserad med hjälp av etylen Etilen Oksitle sterilize edilmiştir
	Non-pyrogenic غير مؤكّد للحمى Непирогенно Neryogenní Ikke-pyrogen Pyrogenfrei Μη πυρετογόνο Apirógeno Mittepürogeenne Pyrogeeniton Non pyrogène Nezapaljivo Nem pirogén Apirogeno Nepirogeninis Niet-pyrogeen Ikke-pyrogen Produkt niepirogenny Não pirogénico Apirogen Апирогенно Neryogénne Apirogeno Nezapaljivo Ikke-pyrogen Pirojenik değildir

Do not reuse | لا يُستخدم ثانية | За еднократна употреба | Nepoužívejte opakovaně | Må ikke genbruges | Nicht wiederverwenden | Να μην επαναχρησιμοποιείται | No reutilizer | Ärge kasutage korduvalt | Ei saa kasutada uudele | Ne pas réutiliser | Nije za ponovnu upotrebu | Nem újrafelhasználható | Non riutilizzare | Nenaudoti pakartotinai | Niet opnieuw gebruiken | Ikke til gjenbruk | Nie używać ponownie | Não reutilizer | A nu se reutiliza | Не использовать повторно | Opakovane nepoužívejte | Samo za enkratno uporabo | Nemojte koristiti više puta | Återanvänd inte | Tekrar kullanmayın

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1 Device Description

SeQure™ Snare System, including General Snare System, Petite Snare System and Micro Snare System according to different loop sizes and the length of the snare, consists of a snare, a snare catheter, a snare loader and a torquer. The snare is constructed of a nitinol shaft and a loop, and they are connected together on the tapered section of the nitinol shaft by the compression force of PTFE heat shrink tubing. The loop is made of braided nitinol wires coiled with a Titanium Nitride coated tungsten wire to enhance radiopacity. The plane of the loop is perpendicular to the nitinol shaft. The pre-formed snare loop can be introduced through catheter without risk of deformation because of shape memory and the super-elastic properties of nitinol. The snare catheter is made of Polyethylene (PE) tubing. The Luer Lock connector has a radiopaque marker band at its distal tip to enhance the visibility during the operation. The snare is pre-assembled with the snare loader and the torquer, and then packaged with corresponding catheter. (See Figure 1 to Figure 5)


Figure 1 SeQure™ Snare System


Figure 2 Snare


Figure 3 Snare catheter


Figure 4 Snare loader


Figure 5 Torquer

2 Indications and Usages

- **Intended users:**
The SeQure™ Snare System is intended to be used in an operating room or catheterization lab of hospitals or clinics by professionally trained clinicians.
- **Intended Use:**
The SeQure™ Snare System is intended for use in the peripheral vascular system and cardiovascular system including the heart and blood vessels of the central circulatory system to retrieve and/or manipulate foreign objects.
- **Indications for Use:**
The SeQure™ Snare System is indicated for use in those clinical conditions that require manipulation, retrieval and/or repositioning of intravascular foreign objects such as indwelling central venous access catheters, inferior vena cava filters, coils and guidewire pull-through procedures within the peripheral vascular system and cardiovascular system including the heart, blood vessels of the central circulatory system.
- **Intended Patient target group:**
After clinical evaluation by doctors, patients who are suitable to use the SeQure™ Snare System for foreign body capture and have no obvious contraindications.

3 Contraindications

- This device is not intended to remove foreign bodies that have become entrapped by tissue growth or endothelialization.
- This device should not be used for fibrin sheath stripping in the presence of atrial

or ventricular septal defects or Persistent Foramen Ovale.

- This device is not intended for the removal of implanted pacing leads.
- Patients who cannot tolerate anticoagulation therapy.
- Patients with active infection.

4 Warnings

- Pulling force applied to catheters during fibrin sheath stripping may damage, stretch, or break indwelling catheters (6 French or smaller in diameter).
- Excessive force applied to retrieve and/or manipulate foreign objects may result in damaging the snare system.
- Contents are supplied sterile, using ethylene oxide (EO). Do not use if sterile barrier is damaged. If damage is found, contact the manufacturer.
- For single patient use only. Do not re-use, re-process or re-sterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Do not use substances that are incompatible with any of the product's components. (See product components section)
- After use, dispose of the product and packaging in accordance with hospital, administrative and/or local government policy.

5 Precautions

- Care should be observed when using this device for removal of a large fibrin sheath in order to minimize risk of pulmonary embolism.

6 Potential Complications

- ✧ Thrombus formation
- ✧ Embolism
- ✧ Air embolism
- ✧ Stroke
- ✧ Myocardial Infarction (depending on device position)
- ✧ Pulmonary Embolism
- ✧ Blood vessel perforation
- ✧ Blood vessel injury
- ✧ Inflammation
- ✧ Anaphylaxis
- ✧ Hematoma at the puncture site
- ✧ Infection
- ✧ Chronic tissue injury
- ✧ Bruising at access site
- ✧ Arrhythmia
- ✧ Low blood pressure

- ✧ Tricuspid valve injury
- ✧ Arterial spasms
- ✧ Device entrapment
- ✧ Difficult to Remove
- ✧ Failed surgery

7 Product Features

- Product Specifications

The SeSure™ Snare System is available with specifications in the following (Table 1).

Table 1 Product specifications of the SeSure™ Snare System

Spec.	Loop Dia. (D1)/mm	Snare Length (L1)/cm	Catheter Size (D2)	Catheter Length (L2)/cm
Snare System-General				
LT-SG-05	5	120	4 Fr	100
LT-SG-10	10	120	4 Fr	100
LT-SG-15	15	120	6 Fr	100
LT-SG-20	20	120	6 Fr	100
LT-SG-25	25	120	6 Fr	100
LT-SG-30	30	120	6 Fr	100
LT-SG-35	35	120	6 Fr	100
Snare System-Petite				
LT-SP-10	10	65	4 Fr	50
LT-SP-25	25	65	6 Fr	50
Snare System-Micro				
LT-SM-02	2	200	3 Fr	175
LT-SM-04	4	200	3 Fr	175
LT-SM-07	7	200	3 Fr	175

Note: 1Fr=1 French=1/3 mm

- Production Components

The following components may come in contact with skin or blood:

SeSure™ Snare System	
Components name	Material
Snare	Nitinol Tungsten with Titanium Nitride coat PTFE
Snare catheter	PE HDPE Tantalum Marker band

SeSure™ Snare System	
Components name	Material
	Barium Sulfate
Loader	PE HDPE Barium Sulfate
Torquer	POM Brass

8 Directions for Use

- **Assembling the Snare**
 1. Select the appropriate size SeSure™ Snare System for the location of the foreign body. In general, the snare loop should be approximately the size of the vessel in which it will be used.
 2. Remove the snare and snare catheter from their protective holders and inspect for any damage.
 3. Remove the snare loader and Torquer from the proximal end of the snare shaft.
 4. Load the snare into the snare catheter by inserting the proximal (non-looped) end of the snare into the distal (non-hubbed) end of the snare catheter, until the proximal end of the snare shaft exits the distal hub and the loop end can be retracted into the distal end of the snare catheter.
 5. Inspect the device by extending and retracting the snare loop through the distal end of the catheter 2-3 times, while carefully examining the catheter, radiopaque marker band and snare for defects or damage.
 6. When appropriate, the snare/snare catheter can be advanced into the vessel as one unit, assembled as described above.
- **Alternative Procedure**
 1. If the snare catheter is already positioned within the vasculature, the provided loader (located on the proximal end of the snare and just distal to the steering handle) may be used to position the snare loop in the indwelling snare catheter. To do so:
 2. Remove the snare from the protective holder and inspect for any damage.
 3. Move the provided loader (located on the proximal end of the snare and just distal to the steering handle) distally until the loop of the snare is enclosed within the tubing portion of the loader.
 4. Insert the distal end of the loader into the hub of the indwelling snare catheter until resistance is felt. This will indicate the tip of the loader is properly aligned with the inner lumen of the catheter.
 5. Keeping the loader as straight as possible, grasp the shaft of the snare just behind the hub of the loader and advance the snare until it is secure within the lumen of the catheter. The loader can be removed at this time by first removing

the steering handle and pulling the loader off the proximal end of the snare's shaft.

- **Snare-Assisted Manipulation/Retrieval**

- 1.Remove the indwelling balloon or delivery catheter, if present. It may be necessary to change or extend the indwelling guidewire to facilitate balloon removal and to upsize the indwelling guide catheter to accommodate the snare.
- 2.If a guidewire is indwelling in the patient at the location of the foreign body, cinch the snare over the proximal end of the wire and advance the snare into the guide catheter or sheath until the distal end of the snare catheter is positioned just proximal to the foreign body.
- 3.If no guidewire is present, then pull the snare just into the distal end of the snare catheter and advance in the guide catheter until it is positioned just proximal to the foreign body.
- 4.Push the snare shaft forward gently to completely open the loop. The loop is then slowly advanced forward, and around the proximal end of the foreign body.
- 5.By advancing the snare catheter, the loop of the snare is closed to grab the foreign body. (Note that attempting to close the loop by pulling the Snare within the Snare catheter will move the loop from its position around the foreign body.)
- 6.To manipulate a foreign body, maintain tension on the snare catheter to retain the hold on the foreign body, and move the snare and snare catheter together to manipulate the foreign body to the desired position.
- 7.To retrieve a foreign body, maintain tension on the snare catheter and move the snare and snare catheter together (proximally) to or into the guide catheter or sheath. The foreign body is then withdrawn through or together with the guiding catheter or vascular sheath. Withdrawal of large foreign bodies may require the insertion of larger sheaths, or guiding catheters, or a cut down at the peripheral site.

- **Snare-Assisted Removal of Fibrin Sheaths from Indwelling Catheters**

- 1.From a common femoral vein approach, advance the selected snare to the inferior vena cava or right atrium.
- 2.Advance a 0.035 inch (0.89mm) guidewire through the port (distal or venous port if more than one lumen) of the indwelling catheter and into the right atrium or inferior vena cava.
- 3.Encircle guidewire with the snare loop.
- 4.Advance the snare over the distal end of the catheter to a position proximal to the fibrin sheath.
- 5.Close the snare around the catheter and continue applying light traction while gently pulling the snare down (towards distal end of catheter) over the port(s).
- 6.Repeat steps 4&5 until the catheter is free of fibrin sheath.

- **Snare-Assisted Venous Canalization**

- 1.Introduce the snare at a patent venous access site and position in the vasculature at the desired venous access site.
- 2.Open snare loop to provide a target to guide a puncture into the desired venous access site.
- 3.Introduce guidewire through needle and through snare loop.
- 4.Close snare over needle by advancing snare catheter.
- 5.Pull the guidewire into inferior vena cava.

9 Clinical Benefit

Clinical benefit is determined by the indication of use of the device as well as clinical background information on the condition treated. A rare but serious complication of endovascular therapy is the unintentional loss or misplacement of an endovascular device within the vasculature. The snares are intended for use in the cardiovascular system as a tool to retrieve and manipulate foreign objects, without having a direct therapeutic or diagnostic function themselves. They are used in a wide range of clinical procedure scenarios. Due to the snares has fracture resistance and superiority of manipulation and retrieval, and minimally invasiveness because of percutaneous, transcatheter method reduce the potential for trauma during punctures while retrieving and manipulating foreign objects.

The technical success could be achieved in most of the patients using snares. Various series have reported a technical success rate of 87% to 100%. (Referenced Literature)

The state of art show that the complications occurrence rate of snares not exceed 6%, though that may be unrelated with device/procedure.

Reference Literature: Wolf, F., Scherthaner, R.E., Dirisamer, A. et al. Endovascular Management of Lost or Misplaced Intravascular Objects: Experiences of 12 Years. Cardiovasc Intervent Radiol 31, 563–568 (2008).

10 Expiry Date

The SeSure™ Snare System is sterilized by ethylene oxide. The shelf life is indicated on the label, and the production date and expiry date are marked on the label. Do not use expired product.

11 Storage and Disposal

- The SeSure™ Snare System should be kept in a protective, well ventilated and clean room. Keep away from direct sunlight, keep dry, store at room temperature.
- Do not squeeze, expose to liquids, bend or expose to extreme heat during transportation.
- If there are any questions before using the SeSure™ Snare System, contact the manufacturer directly.
- After use, dispose of the device and packaging in accordance with hospital, administrative and/or local government policy.

12 Package and Label

The SeQure™ Snare System is supplied sterile. All the components of the SeQure™ Snare System are fixed and protected with protective holders, then sealed in two Tyvek1073B bags, to which a primary label is applied. The product is sterilized and put in a box with IFU and customer feedback form. The primary label is applied to the box. Please contact Lifetech for more information on Instructions for use. (website address: <http://www.lifetechmed.com>, tel: 86 755 86026250).

13 Adverse event feedback

Any complications related to the SeQure™ snare system should be reported to the customer service department of Lifetech Scientific (Shenzhen) Co., Ltd whose mail address is lifetechmed@lifetechmed.com.

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