

# The ALN OATF

Optional Vena Cava Filter

*Just like the original,  
talented and reliable*

**NEW**

**OATF**

## The ALN OATF Optional Vena Cava Filter with and without hook

Based on the expertise  
and experience of ALN

- The only optional IVC Filter in Stainless steel without welding points: 3T MRI Compatible
- 4 approaches
- Approved with an indicated IVC diameter up to 28 mm
- Enhanced stability with its 3 anti-tilting legs

## Optimal Design

3

anti-tilting

6  
legs  
anchoring  
legs

3

centering  
legs

➤ A design with an unique technology:

- Anti-tilting legs for an optimum centering
- Anchoring legs with an anti-transfixion loop to minimize the risk of cava penetration
- Centering legs

➤ The ALN OATF optional IVC Filter with hook can be removed with one of our three extraction kits and/or a lasso

## THE ALN OATF OPTIONAL VENA CAVA FILTER 4 Approches

Femoral  
&  
Popliteal

3T MRI  
compatible

Jugular  
&  
Brachial

## PRODUCT CODES

		Product Code	Maximum diameter of the vena cava (mm)	Length of the filter (mm)	Introducer Sheath	
					French size (ID)	Length (mm)
ALN OATF Vena Cava Filter with Hook	Jugular Kit	FJH.OATF19	28	59	7	600
	Brachial Kit	FBH.OATF19	28	59	7	1200
	Femoral Kit	FFH.OATF19	28	59	7	600
	Popliteal Kit	FPH. OATF19	28	59	7	1200

The implantation kit includes : the ALN filter in its holder, an introduction system (dilator, introducer sheath, pushing catheter), a puncture needle, a «J» Tip guidewire.