

# ALN<sub>2IN1</sub>

## Extraction Kit

Product code : RK-2010

### Unique 2 in 1 RS

ALN Vena Cava Extraction  
and/or Repositioning Kit

Only by jugular approach

Ability  
to change  
the angle  
of the pincer

The choice  
is yours!

A different option to  
remove ALN IVC Filters  
using a snare or the 2IN1  
Extraction Kit

**95.5%**  
Successful  
removal rate  
using the ALN  
Extraction Kit\*

Possibility  
to reposition  
the filter  
if needed


- Approved by the FDA to remove or reposition ALN IVC Filters
- 9F ID Extraction Sheath
- The kit is complete and includes: A pincer catheter, an introduction system, a puncture needle, a "J" guidewire
- Excellent maneuverability
- Safety of the procedure due to the pincer design
- High visibility of the device during the procedure



High



Moderate



Straight

*When the 9F ID Extraction Sheath is lowered by 1 cm, the angle of the curved section decreases. When the 9F ID Extraction Sheath is lowered by 1.5 cm, the ALN 2IN1 RS is straight.*